FUTURO

performance and perfection

152380.0100

HDMI Kamera FUTURO AUTOFOCUS

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HDMI Camera FUTURO AUTOFOCUS

Operating Instructions

CAMÉRA HDMI FUTURO AUTOFOCUS

Mode d'emploi

FOTOCAMERA HDMI FUTURO AUTOFOCUS

Manuale d'uso

CÁMARA HDMI FUTURO AUTOFOCUS

Manual de instrucciones

HDMI KAMERA FUTURO AUTOFOCUS

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KAMERA HDMI FUTURO AUTOFOCUS

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CAMERĂ VIDEO HDMI FUTURO AUTOFOCUS

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HDMI KAMERA FUTURO AUTOFOCUS

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1. Grundlegende Funktionen

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www.brw.ch · info@brw.ch

1. Basic functions

This camera is a CMOS camera with multiple interfaces (HDMI+WLAN+SD card, "X" stands for multiple interfaces) and autofocus function (F stands for autofocus). The image selector of the camera is an ultra-high performance CMOS sensor from Sony. HDMI + WLAN are used as a data transmission interface to an HDMI screen or computer.

When using the HDMI output, XCamView is loaded and a camera control panel and a toolbar are displayed on the HDMI screen. In this case, the USB mouse can be used to set the camera, search and compare captured images and play videos.

The HDMI output can be used to easily obtain a sharp image by means of the integrated automatic / manual focus function. Manual focusing via the coarse / fine adjustment knob on the microscope is not necessary.

To use the WLAN output, disconnect the mouse from the device, plug in the USB WLAN adapter and connect the WLAN port on the computer to the camera. The video stream can now be transferred to the computer using the state-of-the-art software.

Figure 1

1. Basic functions

Basic functions of the camera:

- All-in-one camera (HDMI+WLAN) with C-mount and highly sensitive CMOS sensor from Sony;
- **Automatic / manual focus with sensor movement;**
- For HDMI application, with integrated XCamView software in several languages. The camera function can be controlled with the USB mouse via XCamView. XCamView can also perform other basic processing and control functions;
- 1920 × 1080 (1080p) resolutions to match the current HD displays on the market; support for plug-and-play applications;
- For an HDMI application, an image with a resolution of 2.0 M (1920 x 1080) can be recorded and saved; videos can be recorded and saved as 1080P video stream (asf format);
- With the USB WLAN adapter, the camera can be used as a WLAN camera. The state-of-the-art image processing software is used for playing videos and displaying photos. Support for plug-and-play applications;
- Ultra-fine colour engine with perfect colour reproduction capability (WLAN);
- With advanced video and image processing application software, including professional image processing functions such as 2D measurement, HDR, image stitching, EDF (Extended Depth of Focus), image segmentation and image counter, image stacking, colour composite and noise reduction (USB);
- The camera meets the requirements of various applications and can be used in many areas of industrial inspection, education and research, material analysis, precision measurement, medical analysis, etc.

Possible applications for the camera:

- Scientific research, education (teaching, demonstration and academic exchange);
- Digital laboratories, medical research;
- Industrial visualisation (PCB testing, IC quality control);
- Medical treatments (pathological observation);
- Food products (observation and determination of the total bacteria count in microbial colonies).

1. Basic functions

1.1 Data sheet

Order no.	Sensor and size (mm)	Pixels (μm)	G sensitivity Dark signal	FPS / Resolution	Binning	Exposure
152380.0100	1080p/2M/Snur IMX185 (C) 1 / 1.9" (7.20 x 4.05)	3.75 x 3.75	1120mv with 1/30 s 0.15 mv with 1/30s	60/1920x1080 (HDMI) 25/1920 x 1080 (WLAN)	1 x 1	0.06 ms ~ 918 ms

C = Colour; M = Monochrome

Interface and key functions		
	USB	USB mouse / USB WLAN adapter
	HDMI	HDMI output
	DC12V	12 V/1 A current input
	SD	SD card slot
	ON/OFF	On/off switch
	LED	Power indicator
Additional specification for HDMI output		
UI operation	With USB mouse for use with the embedded XCamView	
Image capture	JPEG format with a resolution of 2 M on SD card	
Video recording	ASF format 1080p 30fps on SD card (8G)	
Camera control panel	Including exposure, gain, white balance, colour correction, sharpness and noise reduction control	
Toolbar	Including zoom, flip, compare, still image, crosshairs, browser function, multiple languages selectable, XCamView version information	
Additional specification for WLAN output		
UI operation	Software on Windows/Linux/OSX/Android platform	
WLAN performance	802.11n 150 Mbps; RF output power 20 dBm (maximum)	
Maximum number of connected devices	3 ~ 6 (depending on environment and connection distance)	
White balance	Automatic white balance	
Colour technology	Ultra-fine TM colour engine (WLAN)	
API for image recording / control	Standard SDK for Windows/Linux/Mac (WLAN)	
Recording system	Still image or film (WLAN)	

Software environment (for USB 2.0 connection)	
Operating system	Microsoft® Windows® XP / Vista /7/8/8.1/10 (32 and 64 bit) / OSx (Mac OS X) / Linux
PC requirements	CPU: Equivalent to Intel Core2 2.8 or higher
	Random access memory: 4 GB or more
	USB port: USB 2.0 high-speed connection (as power connection only, not for data transfer via USB)
	Display: 19" or larger
	CD-ROM
Operating environment	
Operating temperature (in centigrade)	-10~50
Storage temperature (in centigrade)	-20~60
Operating humidity	30~80 % RH
Storage humidity	10 ~ 60% RH
Power supply	12 V DC / 1 A adapter

1.2 Camera and microscope

Figure 2 and rear side

Figure 3 different views

Figure 4 with microscope

1.1.1 Dimensions


Figure 5 Dimensions

1.2 Camera and microscope

Figure 6 Packaging information for

Standard packaging list			
A	Gift box: L: 25.5 cm W: 17.0 cm H: 9.0 cm (1 pc., 1.43 kg/box)		
B	HDMI camera		
C	Power unit: Input: AC 100~240 V 50 Hz/60 Hz, output: 12 V DC / 1 A American Standard: Model: GS12U12-P1I 12W/12 V/1 A: UL/CUL/BSMI/CB/FCC EMI Standard: EN55022, EN61204-3, EN61000-3-2, -3, FCC Part 152 Class B, BSMI CNS14338 EMS Standard: EN61000-4-2, 3, 4, 5, 6, 8, 11, EN61204-3, Class A Light Industry Standard European Standard: Model: GS12E12-P1I 12W/12V/1A; TÜV (GS)/CB/CE/ROHS EMI Standard: EN55022, EN61204-3, EN61000-3-2, -3, FCC Part 152 Class B, BSMI CNS14338 EMS Standard: EN61000-4-2, 3, 4, 5, 6, 8, 11, EN61204-3, Class A Light Industry Standard		
D	HDMI cable		
E	USB mouse		
F	Wireless network adapter with USB interface		
G	CD (driver and utility programs, Ø12 cm)		
Optional accessories			
I	Fixed lens adapter	C-mount for eyepiece tube with 31.75 mm diameter	108005/FMA037
J	108015 (diameter 23.2 mm to 30.0 mm ring) / adapter rings for 30 mm eyepiece tube		

1. Basic functions

1.1.3 Extension of the camera with a microscope or telescope adapter

Extension	Image	
C-mount camera		Image processing; medical imaging; semiconductor equipment; test equipment; document scanner; 2D barcode reader; web camera and security video; microscope imaging;
Microscope camera	 XFCAM1080PHD+FMAXXX(23.2mm Adapter)	

2. Functions on the rear side of the camera


Figure 7 Layout of the rear side of the camera

This is an HDMI camera with autofocus function. By precisely controlling the sensor position, the image can be automatically focused for a stereo or biological microscope. However, this autofocus principle contradicts the principle of image-forming conjugate planes applicable to the microscope and we believe that only a small focus setting is necessary to maintain high image quality.

With a basic online observation, however, the work efficiency can be increased considerably with the camera and manual focusing is no longer required.

3. Quick start guide for the camera

Before switching on the camera, please connect the standard C-mount camera to the camera adapter and connect this to the third microscope tube, which transfers the median image of the microscope object to the camera sensor.

3.1 WLAN mode

1. Plug the **12V/1A** power cable into the **power supply interface** ③ to connect the camera to the power supply. The **LED display** ⑥ then lights up red;
2. Press the **ON/OFF** button ⑤ to switch on the camera. The **LED display** ⑥ then lights up blue;
3. Insert the WLAN antenna included in the scope of delivery into the **WLAN ANTENNA/USB PORT** ② to generate the **WLAN** signal;
4. As soon as the indicator on the **WLAN** antenna begins to flash, connect the computer (or an iPad or smartphone) to the **WLAN** signal, the name of which begins with **BHDC-AFC202M**. The **password** is 12345678;
5. Open the software and start by clicking the name of the camera model in the camera list.

3.2 HDMI mode

1. Connect the HDMI cable to the **HDMI port** ① to connect the camera to the HDMI display;
2. Connect a USB mouse to the **USB port** ② to control the camera with the integrated **XCamView** software;
3. Plug the **12 V/1 A** power unit into the **power supply interface** ③ to connect the camera to the power supply. The **LED display** ⑥ then lights up red;
4. Insert the SD card into the **SD card slot** ④ to save captured images and recorded videos;
5. Press the **ON/OFF** button ⑤ to switch on the camera. The **LED display** ⑥ then lights up blue;
6. Move the mouse pointer to the left side of the video window. The **camera control panel** is then displayed. The control panel provides the functions **manual / automatic exposure, white balance, sharpness, noise reduction**, etc. For details, see 3.3.1
7. When you move the mouse pointer to the top of the video window, a **measurement toolbar** with calibration and other measurement tools is displayed. For further information, see 3.3.3. The measurement data can be exported as a ***.CSV** file.
8. When you move the mouse pointer to the bottom edge of the video window, a **camera control toolbar** is displayed. This enables operations such as **zooming in, zooming out, flipping, still images, crosshairs, WDR** etc. to be performed. For more information, see 3.3.2
9. When you move the mouse pointer to the bottom edge of the video window, the **camera control toolbar** opens automatically. When you click the button **AF**, the **autofocus control panel** appears, which allows you to perform autofocus operations; 3.3.4 and 3.3.5

3.3 Brief introduction to the camera and its functions

The camera UI in Figure 8 shows the **camera control panel** on the left side of the video window, the **measurement toolbar** at the top of the video window, the **camera control toolbar** at the bottom of the video window, and the **autofocus control panel** on the right side of the video window.



Figure 7 Layout of the rear side of the camera

Notes	
1	When the mouse pointer is moved to the left side of the video window, the camera control panel opens automatically;
2	When the mouse pointer is moved to the bottom edge of the video window, the camera control toolbar opens automatically;
3	When the mouse pointer is moved to the bottom edge of the video window, the camera control toolbar opens automatically. Click the button AF to open the autofocus control panel for autofocus operation;
4	Moving the mouse pointer to the top edge of the video window opens the measurement toolbar for the calibration and measurement processes. When the Move / Fixed button is clicked in  of the measurement toolbar , the measurement toolbar remains fixed in one position. In this case, the camera control panel does not open automatically, even if the mouse pointer is moved to the left side of the video window. Only when the button  in the measurement toolbar is clicked to end the measurement process can other operations be carried out in the camera control panel , autofocus control panel or camera control toolbar . When a specific measurement object is selected during the measurement process, an object position and attribute control bar  is displayed, which can be used to change the position and the properties of the selected objects.

3.3.1 The camera control panel on the left side of the video window

Camera control panel	Function	Description of functions
	Snapshot	Record an image or take a snapshot in the currently open video window.
	Video recording	Record a video in the currently open video window.
	Auto Exposure	When automatic exposure is enabled, the system automatically adjusts the exposure time according to the exposure compensation value.
	Exposure compensation	Available when the automatic exposure is enabled. Move left or right to adjust the exposure compensation according to the current video brightness and obtain the appropriate brightness value.
	Exposure time	Available when the automatic exposure is disabled. Move left or right to increase or decrease the exposure time to adjust the brightness of the video.
	Gain	Adjust the gain to increase or decrease the video brightness. The noise is reduced or amplified accordingly.
	Red	Move left or right to increase or decrease the red component in the video window.
	Green	Green is a reference base and cannot be changed.
	Blue	Move left or right to increase or decrease the blue component in the video.
	White balance	Automatic white balance according to the video window.
	Sharpness	Adjust the image sharpness of the video window.
	Noise reduction	Adjust the noise reduction level of the video window
	Saturation	Adjust the saturation level of the video window.
	Gamma	Correction of the gamma level of the video. Move right to increase and left to decrease the gamma level.
	Contrast	Correction of the contrast value of the video. Move right to increase and left to decrease the contrast.
	DC	For DC lighting, there is no fluctuation in the light source, so there is no need to eliminate flicker.
	AC (50 HZ)	Check AC (50 Hz) to eliminate the "flickering" that occurs when lighting at 50 Hz.
	AC (60 HZ)	Check AC (60 Hz) to eliminate the "flickering" that occurs when lighting at 60 Hz.
	Default settings	All settings in the camera control panel are set to the default values.

The camera is controlled via the **camera control panel** to achieve the best image quality for each specific application. The control panel opens automatically when the mouse pointer is moved to the left side of the video window (the **camera control panel** does not open in measurement mode). The **camera control panel** only opens when the measurement process is finished, when the mouse pointer is moved to the left side of the video window). Clicking the button activates the Show / Auto Hide button on the **camera control panel**;

3.3.2 Icons and functions of the camera control toolbar at the bottom edge of the video window

Icon	Function	Icon	Function
	Enlarge the video window		Shrink the video window
	Flip horizontal		Flip vertical
	Video still image		Show crosshairs
	WDR		Start the autofocus control panel
	Search images and videos on the SD card		Settings
	Check version of XCamView		

The setting function  is a little more complex than the other functions. More information about this function is provided below:


Figure 9 Detailed page for WLAN settings

Channel: WLAN signal **channel**. Avoidance of interference caused by using the same channel. Recommendation to choose different channels for different cameras when several WLAN cameras are operated simultaneously;

SSID: Name of the WLAN signal. Can be personalised using the soft keyboard below;

Password: **Password** for the WLAN signal. The **password** can be personalised using the soft keyboard below;

Default settings: Reset **channel**, **SSID** and **password** to the default values;


Figure 10 Comprehensive measurement settings page

Global: Used to set the digits after the decimal point in measurement results;

Calibration line width: Used to specify the line width during measurement and calibration;

Colour: Used to specify the colour of the lines during measurement and calibration;

End point type: Used to specify the shape of the line end points during measurement and calibration: "Null" means no end points, "Rectangle" means rectangular end points. Simplifies calibration;

Point, angle, line, horizontal line, vertical line, rectangle, circle, ellipse, annulus, two circles, polygon, curve: By clicking  next to the measurement patterns listed above, the corresponding attribute settings open to specify the individual properties of the measurement object.


Figure 11 Comprehensive page for units of measurement, calibration and magnification

Name: Names such as 4X, 10X, 20X, 40X, 100X refer to the magnification of the microscope. For microscopes with infinitely variable zoom, it is important to ensure that the selected magnification corresponds to the scale alignment line;

Resolution: Pixels per metre. Devices such as microscopes offer a high resolution;

Delete all: Delete all calibrated magnifications and resolutions;

Delete: Click **Delete** to delete the selected element for a specific resolution;


Figure 12 Page for setting the image format

JPEG: Save the captured image on the SD card in **JPEG** format;

TFT: Save the captured image on the SD card in **TFT** format; in the **TFT** format, not only the image data are stored, but also the measurement data for the image. The camera control and image processing software can open a **TFT** file;


Figure 13 Detailed page for setting the video coding

MJPEG: For saving recorded videos coded in **MJPEG** format;

H264: For saving recorded videos coded in **H264** format;


Figure 14 Detailed page for setting the SD card

Current file system: The maximum storable size of a **FAT32** file is 4 GB; for **NTFS** files this is 2.048 GB. Proposal to convert **FAT32** files to **NTFS** format on a computer; **status unknown:** SD card or file system not recognised;

Figure 15 XFCAM Comprehensive page for specifying the language selection settings

- English:** Set English as the language for all software;
- Chinese (simplified):** Set Chinese (simplified) as the language for all software;
- Chinese (traditional):** Set Chinese (traditional) as the language for all software;
- Korean:** Set Korean as the language for all software;
- Thai:** Set Thai as the language for all software;

3.3.3 Measurement toolbar at the top edge of the video window

The **measurement toolbar** opens when the mouse pointer is moved near the top edge of the video window. The various functions of the **measurement toolbar** are explained below:


Figure 16 Button for measurement toolbar at the top edge of the video window


Icon	Function
	Move / Fix button of the measurement toolbar
☒ Visible	Specify the measurement object in the Show / Hide mode
Pixel	Selection of the desired unit of measurement
NA	To ensure the accuracy of the measurement results when the unit of measurement is not pixels, select the same magnification as the microscope
	Object selection
	Point
	Angle
	Arbitrary line
	Parallel
	Horizontal line
	Vertical line
	Rectangle
	Circle
	Ellipse
	Annulus
	Two circles and centre distance
	Bow
	Polygon

Icon	Function
	Curve
	Perform calibration to determine the appropriate relationship between magnification and resolution, which provides the corresponding relationship between unit of measurement and sensor pixel size. Calibration must be performed using a micrometer. The individual calibration steps can be found in the software manual.
	Conjugated Correction: Click  to execute the Conjugated Correction before performing a calibration. Next, manually adjust the microscope knob for the coarse and fine focus settings, in order to focus the video. Make sure that the magnification in the software matches the magnification of the microscope and then select the appropriate unit of measurement .
	Export the measurement data to a CSV file (*.csv)
	Delete all measurement objects
	Adjustment
	Close the current measurement mode
	Once the measurement is complete, click on a single measurement object to open the object position and property control bar . The icons in the control bar mean: Move left , move right , move up , move down , colour matching and delete .

Note:

- 1) When the user clicks the **Show / Hide** button in  on the **measurement toolbar**, the **measurement toolbar** is fixed. In this case, the **camera control panel** does not open automatically, even if the mouse pointer is moved to the left side of the video window. Only when button  in the **measurement toolbar** is clicked to exit the measurement mode can other operations be carried out in the **camera control panel**, the **autofocus control panel** or the **camera control toolbar**.
- 2) When a specific measurement object is selected during the measurement process, the **object position and attribute control bar**  is displayed, which can be used to change the position and the properties of the selected object.
- 3) To ensure the accuracy of the measurement, click the button **Conjugated Correction**  to reset the camera sensor to the **C-mount** default position before calibration. Measurements can be started after calibration is complete and the video is in focus.
- 4) If the calibration is complete but the camera sensor is not in the **C-mount** position, the **Conjugated Correction** should be performed to reset the sensor to the C-mount default position and the video should be focused before starting the measurement.

3.3.4 Autofocus control panel on the right side of the video window

	Autofocus	When the autofocus button is activated, the system automatically starts the autofocus according to the status of the object until it is in focus;
	Manual focus	When the manual focus function is activated, the position of the camera sensor should be reset by scrolling up or down with the mouse until the object is in focus;
	One push AF	A one-time autofocus process can be performed by clicking the One Push button;
	Conjugated Correction	The camera sensor can be reset to the C-mount default position by clicking the Conjugated Correction button. The Conjugated Correction function allows the sensor position to be calibrated while both the camera video window and the image seen in the eyepiece are in focus. The use of the Conjugated Correction function is recommended when the camera is used for the first time. This will ensure that the camera sensor is in the C-mount default position. This guarantees that the object level, eyepiece image level and camera adapter image level are in the default position; Note: 1) If the height of the object changes, it is important to ensure that the sensor is in the C-mount default position and the focus is set correctly via the coarse and fine adjustment knob of the microscope; 2) The Conjugated Correction should be performed prior to the measurement to guarantee the accuracy of the measurement results (for more information, see Measurement toolbar>Conjugated Correction...).

3.3.5 Focus area in the video window


Figure 17 Focus area

The **focus area** is used to select the area to be examined for the autofocus process. Clicking the button **AF** in the **camera control toolbar** opens the **focus area** together with the **autofocus control panel**. Clicking anywhere in the video window resets the focus area for the **autofocus** process. When the **autofocus panel** is closed, the **focus area** automatically closes.

Note: If the **autofocus** function is active, the **measurement toolbar** will not open when the mouse pointer is moved to the upper video area.

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Brütsch/Rüegger Werkzeuge AG